

Excision of Protein Bands/Spots from Polyacrylamide Gels

1. Purpose

This protocol describes the excision of stained protein spots/bands from a polyacrylamide gel. The excised gel can be submitted for protein identification by mass spectrometry.

2. Reagents and Materials

- 2.1. 1.5 ml (or 0.5 ml) microcentrifuge tubes, siliconized (or low protein-binding)
- 2.2. Acetic acid
- 2.3. Scalpel with fine-sized blades (*e.g.* #11 surgical blade)
- 2.4. Clean work area (*e.g.* positive pressure clean hood or enclosed laminar flow hood)
- 2.5. 90% ethanol
- 2.6. 70% ethanol
- 2.7. Water (Milli-Q or Optima grade, or similar)
- 2.8. Suitable cutting surface (*e.g.* light box, glass plate), very clean
- 2.9. Clean, non-latex gloves

3. Reagent Preparation

- 3.1. 1% Acetic Acid
 - 3.1.1. Add 10 μ l of glacial acetic acid to 990 μ l water (assume 1 ml per 10 gel bands).
 - 3.1.2. Vortex to mix.
- 3.2. Clean tubes
 - 3.2.1. Clean the microcentrifuge tubes with 90% ethanol (or methanol) and dry completely.

4. Procedure

- 4.1. Thoroughly clean the cutting surface with 70% ethanol (or methanol) and place it into the clean work area.
- 4.2. With a clean blade, excise the protein band or region of interest.
 - 4.2.1. If cutting a single band (rather than a multi-band region), do not include any non-stained polyacrylamide. If the band is large, excise only the darkest central portion.
 - 4.2.2. Push the blade straight down to avoid creating rips or tears.
- 4.3. Place the cut gel into a pre-cleaned microcentrifuge tube and add sufficient 1% acetic acid to cover the gel (if commencing an in-gel digest immediately, water can be used instead).
- 4.4. If cutting a large number of spots, add water to the gel as needed to prevent drying. A squirt bottle filled with milli-Q water is handy for this purpose. Dry gel will curl and may rip.
- 4.5. Store at 4°C for short-term or freeze for long-term.

5. Notes

- 5.1. Extreme care must be taken to avoid contaminating the gel with keratin (the major protein of skin and hair) or other environmental proteins. A lab coat with elastic or fitted cuffs should be worn. Wool clothing should not be worn. Clean non-latex gloves must be worn at all times. Prior to cutting, the gel should be kept covered at all times.
- 5.2. Use only the highest-purity reagents for proteomics experiments. Commercially precast gels are often cleaner than lab-made gels.
- 5.3. If the equipment or glassware used to run the gel is also used for Western blots, take great care to avoid cross-contamination from Western blotting reagents, especially “blocking” proteins such as skim milk powder or serum albumin.